

AmScope MD Series Microscope

Quick-Start Guide

Included in the kit:

Microscope with small plastic eyepiece cover
Power cord
USB cord
Microscope lens (10x/18 wide focus) with small rubber cover
Camera (USB attachment for use with a computer) with disc
Green lens for use with light
The World of the Microscope book
10 prepared slides in case
72 blank, glass slides
100 slide covers
1 dropper

1. Plug microscope into a power source using power cable provided, and turn on.
2. Place the microscope lens onto the eyepiece.
3. Place a slide onto the stage and adjust the focus until the object is clear as you look through the lens. This can take some time and patience. Starting with the lowest magnification lens (the red one) is often the best way to start.

For guidance, use “The World of the Microscope” book.

To make your own temporary slides, see page 14.

Using the Camera attachment

1. Remove the microscope lens from the eyepiece, and replace it with the camera piece. Use the USB cord provided to connect the Camera lens to your computer.
2. You will need to use the disc provided or download the software from the Internet.
To download the software, visit
<http://www.amscope.com/software-download>
 - a. Select MD Series as the Camera, and select your operating system from the options provided.
 - b. Choose the most recent version of the software available, and click the “download” link to install it to your computer. Follow the on screen instructions to install.
 - c. Once the install is completed and the software is open, simply select MD35 from the Camera list on the left hand side to see a live image from the microscope.